



| Spare parts | 2000 EVOLUTION | | 2000EVO_SP_01_2006 |
|--------------------|--|------|---------------------------|
| CODE | DESCRIPTION | Q.TY | UNIT PRICE IN EURO |
| 120171/C | MOTOR OPTICAL UNIT including filter disc, motors, photodetector and magnetic hall sensor | 1 | |
| 120172/C | MOTOR UNIT without filter disc & photodetecotor, including magnetic hall sensor | 1 | |
| 120020/C | FILTERS DISC complete with filters | 1 | |
| 120021/C | FILTER DISC with magnet without filters | 1 | |
| 017007 | INTERFERENTIAL FILTER 340 nm | 1 | |
| 017008 | INTERFERENTIAL FILTER 405 nm | 1 | |
| 017009 | INTERFERENTIAL FILTER 505 nm | 1 | |
| 017010 | INTERFERENTIAL FILTER 546 nm | 1 | |
| 017012 | INTERFERENTIAL FILTER 578 nm | 1 | |
| 017011 | INTERFERENTIAL FILTER 630 nm | 1 | |
| 017025 | NEUTRAL DENSITY FILTERS KIT FROM 0.5 TO 2.5 OD (5 pcs) | 1 | |
| 160020/C | PRINTER complete with interface board and wiring | 1 | |
| 150002_1 | FAN 80X80x38 220VAC | 1 | |
| 100070/C | TEMPERATURE SENSOR NTC 5K@25°C | 1 | |
| 132036 | TIP 127 WIRING FOR THERMOSTATE | 1 | |



| Spare parts | 2000 EVOLUTION | | 2000EVO_SP_01_2006 |
|--------------------|---|------|---------------------------|
| CODE | DESCRIPTION | Q.TY | UNIT PRICE IN EURO |
| 080025/C | POLYSNAP POWER ON SWITCH | 1 | |
| 110400/C | MOTHER BOARD | 1 | |
| 080022 | FUSE 500Ma | 2 | |
| 141101/C | TRANSFORMER | 1 | |
| 10040/C | DISPLAY complete with wiring | 1 | |
| 170050/C | LAMP 12V 5W ASSEMBLED | 1 | |
| 110422_2/C1 | PHOTO DETECTOR complete with printed circuit and wiring | 1 | |
| 300101/C | EXTERNAL COVER INCL. KEYBOARD AND DISPLAY | 1 | |
| 300101/V | EXTERNAL COVER | 1 | |
| 100080/C | HALL EFFECT SENSOR complete with printed circuit | 1 | |
| 270001 | MICRO-PROCESSOR MASTER C9 | 1 | |
| 270001 | MICRO-PROCESSOR SLAVE C9 | 1 | |
| 120021/C | FILTER DISC WITH MAGNET | 1 | |
| 017001 | KIT GELATINE (5 pcs for each type) | 1 | |
| 180100 | KEYBOARD FOR 2000EVOLUTION | 1 | |