



Spare parts	Clin Check Plus		CCP_SP_01_2006
CODE	DESCRIPTION	Q.TY	UNIT PRICE IN EURO
120174/C	MOTOR OPTICAL UNIT including filter disc, motors, photodetector and magnetic hall sensor	1	
120172/C	MOTOR UNIT without filter disc and photodetector, including magnetic hall sensor	1	
120020/C1	FILTERS DISC complete with filters	1	
120021/C	FILTER DISC with magnet without filters	1	
017007	INTERFERENTIAL FILTER 340 nm	1	
017008	INTERFERENTIAL FILTER 405 nm	1	
017009	INTERFERENTIAL FILTER 505 nm	1	
017010	INTERFERENTIAL FILTER 546 nm	1	
017012	INTERFERENTIAL FILTER 578 nm	1	
017011	INTERFERENTIAL FILTER 630 nm	1	
017025	NEUTRAL DENSITY FILTERS KIT FROM 0.5 TO 2.5 OD (kit 5 pcs)	1	
100070/C	TEMPERATURE SENSOR NTC 5K@25°C	1	
080025/C	POLY SNAP FEED UNIT POWER ON SWITCH	1	



Spare parts	Clin Check Plus		CCP_SP_01_2006
CODE	DESCRIPTION	Q.TY	UNIT PRICE IN EURO
110202/C	MOTHER BOARD	1	
080022	FUSE 500Ma	2	
141001/C	TRANSFORMER	1	
170040/C	HALOGEN LAMP ASSEMBLED	1	
110422_2/C4	PHOTO DETECTOR complete with printed circuit and wiring	1	
180051	KEYBOARD	1	
100080/C	HALL EFFECT SENSOR complete with printed circuit and wiring	1	
300040_1/V	EXTERNAL COVER FOR CCP	1	
100041/C	DISPLAY complete with wiring for ccp	1	
270005/P	EEPROM PROGRAMMED (2 pcs)	1	
300040_1/C	EXTERNAL COVER FOR CCP WITH DISPLAY AND KEYBOARD	1	
017001	KIT GELATINE (5 pcs for each type)	1	
270001	MICRO-PROCESSOR	1	