



| Biochemical

| Systems

| International Srl

Spare parts	EASY WASH		EW_SP_01_2005
CODE	DESCRIPTION	Q.TY	
150033/C	Compressor	1	
090010	Complete PVC Manifold	1	
120970/S	Well plate Holder	1	
100120/C	Manometer 0-1 bar	1	
400014/C	Air flow regulator	1	
150034	Grid with filter MOD.RCP080V	1	
080036/C	Polysnap	1	
141102	Transformer	1	
400101	Internal tubes set	1	
080031	Fuse 5x20 1.6A rapid	2	
400102	External tubes set	1	
591001/C1	Waste Bottles complete with connectors	1	
591001/C2	Wash Bottles complete with connectors	1	