



Spare Parts	Map Lab & Map Lab Plus		MLP_SP_01_2006
CODE	DESCRIPTION	Q.TY	UNIT PRICE IN EURO
120170/C	MOTOR OPTICAL UNIT including filter disc, motors, photodetector and magnetic sensor hall	1	
120172/C	MOTOR UNIT without filter disc and photodetector, including sensor hall	1	
120020/C3	FILTERS DISC complete with filters	1	
120021/C	FILTER DISC with magnet without filters	1	
017007	INTERFERENTIAL FILTER 340 nm	1	
017008	INTERFERENTIAL FILTER 405 nm	1	
017009	INTERFERENTIAL FILTER 505 nm	1	
017010	INTERFERENTIAL FILTER 546 nm	1	
017012	INTERFERENTIAL FILTER 578 nm	1	
017011	INTERFERENTIAL FILTER 630 nm	1	
017013	INTERFERENTIAL FILTER 492nm	1	
017014	INTERFERENTIAL FILTER 450nm	1	
120857/C	MOTOR OPTICAL UNIT including filter disc for EIA, motors photodetector and magnetic hall sensor	1	
120857/C1	MOTOR OPTICAL UNIT without filter disc for EIA and photodetector, including magnetic hall sensor	1	



Spare Parts	Map Lab & Map Lab Plus		MLP_SP_01_2006
CODE	DESCRIPTION	Q.TY	UNIT PRICE IN EURO
017025	NEUTRAL DENSITY FILTERS KIT FROM 0.5 TO 2.5 OD (kit 5 pcs)	1	
120861/C	FILTERS DISC complete with filters for EIA	1	
080033	FUSE 5X20 1A	2	
110560_4/K	MOTHER BOARD	1	
100070/C	TEMPERATURE SENSOR NTC 5K@25°C	1	
100080/C	HALL EFFECT SENSOR complete with pcb interface	1	
160010/C1	PRINTER complete with printed circuit and wiring	1	
160010	THERMAL PRINTER LTP256	1	
110411_1/K	THERMAL PRINTER INTERFACE	1	
170060/C	LAMP 3.5 V ASSEMBLED for EIA	1	
170020/C	LAMP 12V 20W ASSEMBLED	1	
180185_2	KEYBOARD FOR MAP LAB PLUS	1	
090002_1	PROTECTION COVER FOR CUVETTES	1	
400012_1/C	FLOW CELL 80µL	1	
400030/C	FLOW CELL 18µL with stainless steel join	1	



Spare Parts	Map Lab & Map Lab Plus		MLP_SP_01_2006
CODE	DESCRIPTION	Q.TY	UNIT PRICE IN EURO
150030/C	PERISTALTIC PUMP MOTOR complete with wiring	1	
100040/C2	DISPLAY FOR MAP LAB & MAP LAB PLUS complete with wiring	1	
120210/C	PERISTALTIC ROTOR ASSEMBLED complete with motor	1	
110422_2/C2	PHOTODETECTOR complete with printed circuits and wiring	1	
110422_2/C3	PHOTODETECTOR FOR EIA complete with printed circuits and wiring	1	
141190/C	TRANSFORMER FOR MAP LAB & MAP LAB PLUS complete with wiring	1	
400000	PERISTALTIC PUMP TUBE	1	
400080_1	ASPIRATION ASSEMBLED TUBE FOR FLOW CELL (IN) for stainless steel join	1	
400081_1	ASPIRATION ASSEMBLED TUBE FOR FLOW CELL (IN) with plastic join	1	
400091	CIRCUIT TUBE FLOW CELL/PERISTALTIC TUBE (OUT) for stainless steel join	1	
400091_1	CIRCUIT TUBE FLOW CELL/PERISTALTIC TUBE (OUT) with plastic join	1	
211110	PREAMPLIFIER SIMULATOR BOX (TOOLS)	1	
211100	TOOL BOX	1	
150032/C1	STRIP MOTOR FOR EIA	1	
300032_13	EXTERNAL COVER FOR MLP	1	



Spare Parts	Map Lab & Map Lab Plus		MLP_SP_01_2006
CODE	DESCRIPTION	Q.TY	UNIT PRICE IN EURO
150006	FAN 40X40X20 12VDC	1	
150002_1	FAN 80X80X38 220VAC	1	
017001	KIT GELATINE (5 pcs for each type)	1	
080035	POLY SNAP POWER ON SWITCH	1	