



Spare Parts	Screen Master 3000		SM3000_SP_01_2006
CODE	DESCRIPTION	Q.TY	UNIT PRICE IN EURO
120170/C	MOTOR OPTICAL UNIT including filter disc, motors, photodetector and magnetic hall sensor	1	
120172/C	MOTOR UNIT without filter disc and photodetector, including magnetic hall sensor	1	
120021/C	FILTER DISC with magnet without filters	1	
120020/C	FILTERS DISC complete with filters	1	
017007	INTERFERENTIAL FILTER 340 nm	1	
017008	INTERFERENTIAL FILTER 405 nm	1	
017009	INTERFERENTIAL FILTER 505 nm	1	
017010	INTERFERENTIAL FILTER 546 nm	1	
017012	INTERFERENTIAL FILTER 578 nm	1	
017011	INTERFERENTIAL FILTER 630 nm	1	
017025	NEUTRAL DENSITY FILTERS KIT FROM 0.5 TO 2.5 OD (kit 5 pcs)	1	
160010/C	PRINTER complete with interface board	1	
160010	THERMAL PRINTER LTP256	1	
110411_1/K	THERMAL PRINTER INTERFACE	1	
090001	PROTECTION COVER FOR PRINTER	1	



Spare Parts	Screen Master 3000		SM3000_SP_01_2006
CODE	DESCRIPTION	Q.TY	UNIT PRICE IN EURO
150002_1	FAN 80X80x38 220V AC	1	
100090/C	TEMPERATURE SENSOR LM35CAZ	1	
100000	PELTIER CELL (1 pc)	1	
100000/C	PELTIER CELL ASSEMBLED (2 pcs)	1	
170020/C	LAMP 12V 20W ASSEMBLED	1	
080035	POLYSNAP POWER ON SWITCH	1	
080034	SWITCH 110,220V	1	
110461_4/K	MOTHER BOARD	1	
080033	FUSE 5X20 1A	2	
141121/C	TRANSFORMER	1	
100040/C1	DISPLAY complete with wiring	1	
110422_2/C2	PHOTO DETECTOR complete with printed circuits and wiring	1	
120210/C	PERISTALTIC ROTOR ASSEMBLED including motor	1	
100080/C	HALL EFFECT SENSOR complete with printed circuit	1	
180101	KEYBOARD	1	
270001	MICRO- PROCESSOR MASTER C9	1	



Spare Parts	Screen Master 3000		SM3000_SP_01_2006
CODE	DESCRIPTION	Q.TY	UNIT PRICE IN EURO
270027	MICRO- PROCESSOR SLAVE DALLAS	1	
150030/C	PERISTALTIC PUMP MOTOR complete with wiring	1	
400000	PERISTALTIC PUMP TUBE	1	
090002_1	PROTECTION COVER FOR CUVETTES	1	
400012_1/C	FLOW CELL 80µL	1	
400030	FLOW CELL 18uL for plastic join	1	
400030/C	FLOW CELL 18µL with stainless steel join	1	
100061/C	BATTERIES PACK (3 elements)	1	
400080_1	ASPIRATION ASSEMBLED TUBE FOR FLOW CELL (IN) for stainless steel join	1	
400081_1	ASPIRATION ASSEMBLED TUBE FOR FLOW CELL (IN) with plastic join	1	
400090	CIRCUIT TUBE FLOW CELL/PERISTALTIC PUMP (OUT) for stainless steel join	1	
400090_1	CIRCUIT TUBE FLOW CELL/PERISTALTIC PUMP (OUT) with plastic join	1	
400001	TUBE PVC CRISTAL LIGHT 4x7 (1MT)	1	
300030_6	EXTERNAL COVER FOR SCR3000	1	
150006	FAN 40x40x20 12VDC	1	
017001	KIT GELATINE (5 pcs for each type)	1	
211100	TOOL BOX	1	