



Spare parts	Screen Master Plus		SMP_SP_01_2006
CODE	DESCRIPTION	Q.TY	UNIT PRICE IN EURO
120171/C	MOTOR OPTICAL UNIT including filter disc, motors, photodetector and magnetic hall sensor	1	
120172/C	MOTOR UNIT without filter disc and photodetecotor, including magnetic hall sensor	1	
120020/C4	FILTERS DISC complete with filters	1	
120021/C	FILTER DISC with magnet without filters	1	
017007	INTERFERENTIAL FILTER 340 nm	1	
017008	INTERFERENTIAL FILTER 405 nm	1	
017009	INTERFERENTIAL FILTER 505 nm	1	
017010	INTERFERENTIAL FILTER 546 nm	1	
017012	INTERFERENTIAL FILTER 578 nm	1	
017011	INTERFERENTIAL FILTER 630 nm	1	
017025	NEUTRAL DENSITY FILTERS KIT FROM 0.5 TO 2.5 OD (5 pcs)	1	
160020/C	PRINTER complete with interface board and wiring	1	
150002_1	FAN 80X80x38 220VAC	1	



Spare parts	Screen Master Plus		SMP_SP_01_2006
CODE	DESCRIPTION	Q.TY	UNIT PRICE IN EURO
100070/C	TEMPERATURE SENSOR NTC 5K@25°C	1	
132036	TIP 127 WIRING FOR THERMOSTATE	1	
080025/C	POLYSNAP POWER ON SWITCH	1	
110400/C	MOTHER BOARD	1	
080022	FUSE 500Ma	2	
141101/C	TRANSFORMER	1	
100020/C	DISPLAY complete with wiring	1	
170050/C	LAMP 12V 5W ASSEMBLED	1	
110422_2/C1	PHOTO DETECTOR complete with printed circuit and wiring	1	
300100/CVB	EXTERNAL COVER INCL. KEYBOARD AND DISPLAY (blu)	1	
300100/CVBI	EXTERNAL COVER INCL. KEYBOARD AND DISPLAY (white)	1	
300100/VBI	EXTERNAL COVER INCL. White	1	
300100/VB	EXTERNAL COVER INCL. BLUE	1	
100080/C	HALL EFFECT SENSOR complete with printed circuit	1	



Spare parts	Screen Master Plus		SMP_SP_01_2006
CODE	DESCRIPTION	Q.TY	UNIT PRICE IN EURO
270001	MICRO-PROCESSOR MASTER C9	1	
270001	MICRO-PROCESSOR SLAVE C9	1	
017001	KIT GELATINE (5 pcs for each tupe)		
180020	KEYBOARD FOR SMP AND VET (blue)	1	
180010	KEYBOARD FOR SMP (black)	1	